Mapping of Merozoite Surface Protein 1-19 Binding Regions in Band 3 Peptide of Red Blood Cell Membrane

Heather Mentzer¹, Lixiao Zeng², Christos G.Takoudis², Athar H. Chishti³, University of South Carolina¹, Chemical Engineering Department², Pharmacology Department³, University of Illinois at Chicago

UIC NSF-REU Program

Acknowledgements

NSF EEC-0453432 Grant and NSF-REU Program 2005

Professor Christos G. Takoudis, Professor Andreas Linninger, Professor Kenneth Brezinsky, Professor Dino Megaridis, Professor Patrick Rousche and Professor David Schneeweis

Professor Athar H. Chishti and Lixiao Zeng Members in Prof. Chishtis' lab and AMReL lab

Outline

- Overall project description
- Information on Malaria
- Description of Merozoite surface protein 1-19 in plasmodium falciparum parasite
- Description of the Band 3 peptide in red blood cell (RBC) membranes
- Biacore optical sensor
 - Surface plasmon resonance
- Results
- Conclusions

Overall Purpose and Project

- Understand the mechanism of Band 3 as a crucial invasion receptor for the malaria parasite.
- From this research new vaccines can be discovered
- Project: Evaluate and analyze the binding sites for *MSP1-19in Band 3 peptide region.
- Experiments: concentration dependence
 - *GST-*5ABC, GST-6A, GST-5, GST-6, GST-5ABC6A, GST-6A5ABC with *Trx-19, *Trx-19A, and Trx-19B

^{*}MSP1-19 is Merozoite Surface Protein 1, 19 is the molecular weight in kDa


^{*}glutathione S-transferase (GST), thioredoxin (Trx)- these are fusion components

^{*5}ABC, 6A, 5, 6, 5ABC6A, 6A5ABC are segments in band 3 peptide.

^{*}Trx-19, Trx-19A, Trx-19B are the MSP1-19, MSP1-19A, MSP1-19B protein

Malaria Statistics

Geographic Distribution of Malaria



- Figure depicts the distribution of malaria. This doesn't mean that areas not designated doesn't get the disease.
- 1 million deaths are attributed to malaria yearly
- Over 40% of the world's population is affected by this disease
- Over 300 million cases that report severe symptoms yearly
- Plasmodium falciparum
 - One of the four strands of malaria
 - This is the most widespread and dangerous of the strands

Center of Disease Control & Prevention http://www.cdc.gov/page.do


Mechanism of the invasion of p. falciparum

- *P. falciparum invasion is possible by the erythrocyte inside the RBC
- Two regions in the band 3 peptide may be the host receptors in the RBC membrane
 - 5ABC, 6A
- P. falciparum Merozoite Surface Protein 1 (MSP1) is a major membrane protein
 - Attaches to the merozoite surface from its C-terminus**
 - MSP1 is essential in the invasion of p. falciparum
 - We use MSP1-19 to see if the protein plays a role in the invasion of RBCs.

**C-terminus is the carboxyl group at the end of the MSP1-19 strand

^{*}Vikas K. Goel; Xuerong Li..al ect. Band 3 is a host receptor binding merozoite surface protein 1 during the Plasmodium falciparum invasion of erythrocytes. PNAS, vol. 100, issue 9, Pg 5164-5169. Aug 2003

Mechanism of the invasion of p. falciparum





^{*} Fairhurst, Rick M.; Baruch, Dror I; al ect. *Abnormal display of PfEMP-1 on erythrocytes carrying haemoglobin C may protect against malaria.* Nature, vol 435, issue7045, Pg 1117- 1121. Jun 2005.

MSP1-19 Diagram 2 4 5 7 8 10 11 12 13 15 16 17 30 kDa 83 kDa 38 kDa 42 kDa N **MSP1**₄₂ MSP1₃₈ **Bacterial** Conserved Region **Expression** (1262-1639)(902-1261)Susceptible Binding sites MSP1₁₉ Not considered (1526-1639)# -the sequence of the MSP1 strand MSP1_{19A} MSP1_{10R}

- OMSP1-19 is used to see whether this protein plays a role in the invasion of RBCs.
- OTruncated MSP1-19A, MSP1-19B will be used in the mapping of the binding sites of MSP1-19

Band 3 Peptide diagram





- •Band 3 peptides are more susceptible to invasion by plasmodium falciparum
- •5ABC is a major segment in the band 3 peptide
- •Test 6A to see how it contributes to the binding

Fujinaga, Jocelyne; al. ect. *Topology of the Membrane Domain of Human Erythrocyte Anion Exchange Protein, AE1.* THE JOURNAL OF BIOLOGICAL CHEMISTRY. Vol. 274, No. 10, pp. 6626–6633, 1999

Biacore Optical Sensor



*Biacore 1000-an early model

- BIA- biospecific interaction analysis
 - Utilizes surface plasmon resonance to perform protein interaction analysis
- *Biosensor technology
 - Information that can be obtained, specificity, affinity, concentration and kinetic parameters
- *Sensor chip:
 - CM5
 - Carboxy methyl dextran
 - Consists of a glass support with a thin layer of gold


^{*}giga.sct.ub.es:800/ ImatgesHQ/biacore.jpg


^{*}www.biacore.com

^{*}Jonsson, U.; Fagerstam, L.; al ect. Real-Time Biospecfic Interaction Analysis using Surface Plasmon Resonance and a Sensor Chip Technology. BioTechniques, Vol 11, issue 5. Nov. 1991 pg. 620-627

Surface Plasmon Resonance

- *Polarized light hits the layer of glass between two indices on the sensor chip
 - Higher- glass surface
 - Lower-running buffer in microchannel.
- Light reflects from higher indices and evanescent wave is formed *
 - Penetrates the lower index
- An intensity dip appears-which is the SPR angle.
 - Angle can change as the response changes
 - Response changes due to the accumulation of mass.






^{*}Biacore System Manual

Results for GST-5ABC and GST6A

Binding of GST5ABC-Trx-19

Experiments:

- GST-5ABC & GST-6A with Trx-19, 19A, 19B
- Strong Bindings seemed to occur
- K_D value of Trx-19A is slightly lower than Trx-19B
- 19A & 19B help contribute to the overall binding of Trx-19.





Parameters for Binding

Bindings Kinetic Constants	KA (1/M)	KD (M)	σ² Chi2 value
GST-5ABC-Trx-19	6.45E+06	1.55E-07	1.67
GST-5ABC-Trx-19A	1.31E+07	7.64E-08	6.07
GST-5ABC-Trx-19B	8.52E+06	1.17E-07	1.77
GST-6A-Trx-19	8.01E+06	1.25E-07	9.34
GST-6A-Trx-19A	4.08E+06	2.45E-07	2.17
GST-6A-Trx-19B	1.24E+07	8.07E-08	7.81

Results for GST-5ABC6A and GST-6A5ABC

- Chimeric proteins- composed of parts
- When chimeric proteins-were used as ligand and Trx19 was the analyte, no binding occurred
- When Trx19A was used as ligand and chimeric proteins were the analyte, strong binding occurred
- Kinetic values for the experiments are close to values from the GST-5ABC experiments





GST-6A5ABC with Trx-19A

Parameters for Binding

Bindings Kinetic Constants	KA (1/M)	KD (M)	σ² Chi2 value
GST-5ABC-Trx-19	6.45E+06	1.55E-07	1.67
GST-5ABC-Trx-19A	1.31E+07	7.64E-08	6.07
GST-5ABC-Trx-19B	8.52E+06	1.17E-07	1.77
GST-6A-Trx-19	8.01E+06	1.25E-07	9.34
GST-6A-Trx-19A	4.08E+06	2.45E-07	2.17
GST-6A-Trx-19B	1.24E+07	8.07E-08	7.81
GST-5ABC6A-Trx-19	1.86E+04	5.39E-05	12.0
GST-6A-5ABC-Trx-19	0.828	1.21	12.3
GST-6A5ABC-Trx-19A	9.49E+07	1.05E-08	2.88

Conclusions

- Binding strengths between the interaction of band 3 peptide and MSP1-19 were found.
- 19A, 19B contributes to the overall binding of Trx-19. MSP1-19 has important binding sites for Band 3 peptide
- GST-6A has similar affinity binding as GST5ABC and the kinetic constants show the 6A is also an important receptor
- Combined 5ABC and 6A region works as receptors in the invasion of P. falciparum